


Bioingentech Genomic RNA Purification

Rev. Date: February 2017


RNA pathogens detection



RNA Extraction from Shrimp

Samples: Hepatopancreas, Gills and Pleopods.

- 1** Cut the samples and fill the tube until approximately 200 uL mark in a 2000 uL tube.
- 2** Macerate the sample with scissor.
- 3** Add 1 mL of CLS to wash the sample.
- 4** Centrifuge to 14.000 RPM for 3 minutes then discard supernatant
- 5** Add 250 µL RNA lysis buffer, homogenize and incubate for 10 min at room temperature.
- 6** Add 50 µL of chloroform.
- 7** Shake tube vigorously by hand for 15 seconds and vortex the sample briefly.
- 8** Incubate for 10 min at -20°C.
- 9** Centrifuge the sample at 12000 x g for 15 minutes.
- 10** Remove 45 µL of upper aqueous phase into a new tube and add 45 µL of Isopropanol, homogenize and incubate for 10 minutes -20°C.
- 11** Centrifuge at 12000 x g for 10 minutes.
- 12** Remove the supernatant (70µL) and wash the pellet with 45 µL of RNA Wash solution
- 13** Centrifuge the tube at 12000 x g for 5 min and discard the wash.
- 14** Resuspend the RNA pellet in RNA rehydration solution (5-10µL)
- 15** Add 2 µL of RNA of sample in (4,5 µL Vet PCR RT PCR Premixture+ 4 µL DNase/RNase free water) *



*See Components PCR Kit: Vet PCR operation Manual Bioingentech, TABLE 1. www.kitpcr.com


RNA Lysis Buffer Bioingentech REAGENT

Reagent is ready to use

A. Preparing whole blood with EDTAK3 samples: RNA isolation procedure

A.1 RNA extracted from 62 μ L of blood, using 250 RNA lysis buffer Bioingentech

- 1** Add 250 μ L RNA lysis buffer per 62 μ L of blood, homogenize and incubate for 10 min at room temperature.
- 2** Add 50 μ L of chloroform.
- 3** Shake tube vigorously by hand for 15 seconds and vortex the sample briefly.
- 4** Incubate for 10 min at -20°C .
- 5** Centrifuge the sample at 12000 x g for 15 minutes.
- 6** Remove 45 μ L of upper aqueous phase into a new tube and add 45 μ L of Isopropanol, homogenize and incubate for 10 minutes -20°C .
- 7** Centrifuge at 12000 x g for 10 minutes.
- 8** Remove the supernatant (70 μ L) and wash the pellet with 45 μ L of RNA Wash solution.
- 9** Centrifuge the tube at 12000 x g for 5 min and discard the wash.
- 10** Resuspend the RNA pellet in RNA rehydration solution (3 μ L).
- 11** Add 2 μ L of RNA of sample in (4,5 μ L Vet PCR RT PCR Premixture+ 4 μ L DNase/ RNase free water) *



*See Components PCR Kit: Vet PCR operation Manual Bioingentech, TABLE 1. www.kitpcr.com


RNA Lysis Buffer Bioingentech REAGENT

Reagent is ready to use

A. Preparing whole blood with EDTAK3 samples: RNA isolation procedure

A.2 RNA extracted from 125 μ L of blood, using 500 μ L RNA lysis buffer Bioingentech

- 1** Add 500 μ L RNA lysis buffer per 125 μ L of blood, homogenize and incubate for 10 min at room temperature.
- 2** Add 100 μ L of chloroform.
- 3** Shake tube vigorously by hand for 15 seconds and vortex the sample briefly.
- 4** Incubate for 10 min at -20°C .
- 5** Centrifuge the sample at 12000 x g for 15 minutes.
- 6** Remove 90 μ L of upper aqueous phase into a new tube and add 90 μ L of Isopropanol, homogenize and incubate for 10 minutes -20°C .
- 7** Centrifuge at 12000 x g for 10 minutes.
- 8** Remove the supernatant (160 μ L) and wash the pellet with 90 μ L of RNA Wash solution.
- 9** Centrifuge the tube at 12000 x g for 5 min and discard the wash.
- 10** Resuspend the RNA pellet in RNA rehydration solution (5 μ L).
- 11** Add 2 μ L of RNA of sample in (4,5 μ L Vet PCR RT PCR Premixture+ 4 μ L DNase/ RNase free water) *



*See Components PCR Kit: Vet PCR operation Manual Bioingentech, TABLE 1. www.kitpcr.com


RNA Lysis Buffer Bioingentech REAGENT

Reagent is ready to use

A. Preparing whole blood with EDTAK3 samples: RNA isolation procedure

A.3 RNA extracted from 250 μ L of blood, using 1 mL RNA lysis buffer Bioingentech

- 1** Add 1 mL RNA lysis buffer Bioingentech per 250 μ L of blood, homogenize and incubate for 10 min at room temperature.
- 2** Add 200 μ L of chloroform.
- 3** Shake tube vigorously by hand for 15 seconds and vortex the sample briefly.
- 4** Incubate for 10 min at -20°C .
- 5** Centrifuge the sample at 12000 x g for 15 minutes.
- 6** Remove 180 μ L of upper aqueous phase into a new tube and add 180 μ L of Isopropanol, homogenize and incubate for 10 minutes -20°C .
- 7** Centrifuge at 12000 x g for 10 minutes.
- 8** Remove the supernatant (340 μ L) and wash the pellet with 180 μ L of RNA Wash solution.
- 9** Centrifuge the tube at 12000 x g for 5 min and discard the wash.
- 10** Resuspend the RNA pellet in RNA rehydration solution (10 μ L).
- 11** Add 2 μ L of RNA of sample in (4,5 μ L Vet PCR RT PCR Premixture+ 4 μ L DNase /RNase free water) *



*See Components PCR Kit: Vet PCR operation Manual Bioingentech, TABLE 1. www.kitpcr.com


RNA Lysis Buffer Bioingentech REAGENT

Reagent is ready to use

B.-Preparing white blood cells: RNA isolation procedure

RNA extracted from white blood cells, using 300µL RNA lysis buffer Bioingentech

- 1** Add 1 mL cell Lysis Bioingentech per 300 µl of blood, homogenize and Centrifuge the sample at 12.000 x g for 5 minutes. Remove the supernatant and then repeat the first step.
- 2** Add 300 µl RNA lysis buffer Bioingentech per pellet, homogenize and incubate for 10 min at room temperature. for 10 min at room temperature.
- 3** Add 100 µL of chloroform.
- 4** Shake tube vigorously by hand for 15 seconds and vortex the sample briefly.
- 5** Incubate for 10 min at -20°C.
- 6** Centrifuge the sample at 12000 x g for 15 minutes.
- 7** Remove 80 µL of upper aqueous phase into a new tube and add 120 µL of Isopropanol, homogenize and incubate for 10 minutes -20°C.
- 8** Centrifuge at 12000 x g for 10 minutes.
- 9** Remove the supernatant (180 µL) and wash the pellet with 150 µL of RNA Wash solution.
- 10** Centrifuge the tube at 12000 x g for 5 min and discard the wash.
- 11** Resuspend the RNA pellet in RNA rehydration solution (10µL).
- 12** Add 2 µL of RNA of sample in (4,5 µL Vet PCR RT PCR Premixture+ 4 µL DNase/RNase free water) *



*See Components PCR Kit: Vet PCR operation Manual Bioingentech, TABLE 1. www.kitpcr.com


RNA Lysis Buffer Bioingentech REAGENT

Reagent is ready to use

C.-Preparing serum samples: RNA isolation procedure


RNA extracted from SERUM, using 300µL RNA lysis buffer Bioingentech

- 1** Centrifuge 1,5 mL serum sample at 12.000 x g for 5 minutes. Remove the supernatant.
- 2** Add 300 µl RNA lysis buffer Bioingentech per pellet, homogenize and incubate for 10 min at room temperature.
- 3** Add 100 µL of chloroform.
- 4** Shake tube vigorously by hand for 15 seconds and vortex the sample briefly.
- 5** Incubate for 10 min at -20°C.
- 6** Centrifuge the sample at 12000 x g for 15 minutes.
- 7** Remove 80 µL of upper aqueous phase into a new tube and add 120 µL of Isopropanol, homogenize and incubate for 10 minutes -20°C.
- 8** Centrifuge at 12000 x g for 10 minutes.
- 9** Remove the supernatant (180 µL) and wash the pellet with 150 µL of RNA Wash solution.
- 10** Centrifuge the tube at 12000 x g for 5 min and discard the wash.
- 11** Resuspend the RNA pellet in RNA rehydration solution (6µL).
- 12** Add 2 µL of RNA of sample in (4,5 µL Vet PCR RT PCR Premixture+ 4 µL DNase/RNase free water) *



*See Components PCR Kit: Vet PCR operation Manual Bioingentech, TABLE 1. www.kitpcr.com

Conditions used and suggested of RT- PCR to obtain complementary DNA (cDNA) RT-PCR-1



Add 1,0 µL of Brig RT-PCR Solution and 1,0 µL of Biotech transcriptase solution RT-PCR-2

