

PRODUCTS USE LIMITATIONS WARRANTY, DISCLAIMER.

Bioingentech products contain chemicals which may be harmful if misused. Due care should be exercised with all Bioingentech products to prevent direct human contact.

Each Bioingentech products is shipped with documentation stating specifications and other technical information. Bioingentech products are warranted to meet or exceed the stated specifications.

Bioingentech's sole obligation and the customer's sole remedy is limited to replacement of products free of charge in the event products fail to perform as warranted. Bioingentech makes no other warranty of any kind whatsoever, and SPECIFICALLY DISCLAIMS AND EXCLUDES ALL OTHER WARRANTIES OF ANY KIND OR NATURE WHATSOEVER, DIRECTLY OR INDIRECTLY, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION, AS TO THE SUITABILITY, PRODUCTIVITY, DURABILITY, FITNESS FOR A PARTICULAR PURPOSE OR USE, MERCHANTABILITY, CONDITION OR ANY OTHER MATTER WITH RESPECT TO BIOINGENTECH PRODUCT. In no event shall Bioingentech be liable for claims for any other damages, whether direct, incidental, foreseeable, consequential, or special (including but not limited to loss of use, revenue or profit), whether based upon warranty, contract, tort (including negligence), or strict liability arising in connection with the sale or the failure of Bioingentech products to perform in accordance with the stated specifications.

Product Information

Kit Name: Bioingentech Genomic DNA Purification Kit

Storage Recommendations: 22-25°C

The Bioingentech Genomic DNA Purification Kit provides a simple, reliable, and rapid method for isolating high-quality total DNA from a wide variety of samples including animal and plant cells and tissue, bacteria and yeast.

The result is high yield, high quality DNA that is suitable for downstream applications such as amplification and restriction digest without further manipulation.

The Kit consist of the following components for purifying total DNA from a wide variety of samples.

1. Cell Lysis Solution
2. Nuclei Lysis Solution
3. Proteinase K
4. RNase Solution
5. Wash Solution
6. Proteinase Precipitation Solution
7. DNA Rehydration Solution

Quality Control

Each lot of Bioingentech Genomic DNA Purification Kit was tested successfully, and the results are consistent lot to lot.

Stability

The Bioingentech Genomic DNA Purification Kit is stable up to the expiration date indicated on the kit label.

The Bioingentech Total RNADetection Kit is a registered trademark of Bioingentech Corporation